


EAR CLEANER

INSTRUCTION



- 1) Open the cover to the battery compartment and install the battery properly.
- 2) Push the silicone nozzle on to the nozzle of the ear cleaner.
- 3) Push the switch to the position of "ON", then put the nozzle into the ear.

! ATTENTION

1. CLEAN THE NOZZLE AND FILTER AFTER USING IT EVERY TIME.
2. CHILDREN SHOULD USE IT UNDER ADULT SUPERVISION.
3. USE 2"AA" 1.5V BATTERIES(NOT INCLUDED).

